

## SYBR Safe™ DNA Gel Stain

The safer ethidium bromide alternative

**S-33100** SYBR Safe™ DNA gel stain, 1 L

**S-33101** SYBR Safe™ DNA gel stain, 4 L

**S-33110** SYBR Safe™ DNA Gel Stain Starter Kit, with 1 L of SYBR Safe™ DNA gel stain (S-33100) and one photographic filter (S-37100)

### Quick Facts

#### Storage upon receipt:

- Room temperature

**Ex/Em:** 280, 502/530 nm, bound to DNA

### Introduction

SYBR Safe™ DNA gel stain has been specifically developed for reduced mutagenicity, making it safer than ethidium bromide for staining DNA in agarose or acrylamide gels. SYBR Safe stain comes as a ready-to-use solution that can be used just like an ethidium bromide solution, and the detection sensitivity with SYBR Safe stain is better than with ethidium bromide. DNA bands stained with SYBR Safe DNA gel stain can be detected using a standard UV transilluminator, a visible-light transilluminator or a laser-based scanner. The stain is also suitable for staining RNA in gels. Bound to nucleic acids, SYBR Safe stain has fluorescence excitation maxima at 280 and 502 nm, and an emission maximum at 530 nm (Figure 1).


#### Benefits

- Increased safety — Tests negative in three mammalian cell-based assays for genotoxicity (Table 1); less mutagenic than ethidium bromide in Ames tests (Figure 2); not classified as hazardous waste under U.S. Federal regulations (Table 2)
- Better performance — Twice as sensitive as ethidium bromide
- Convenient — Ready-to-use in 0.5X TBE
- Fast — No destaining required

### Materials

#### Contents

SYBR Safe DNA gel stain is supplied ready-to-use in two sizes. The 1 L unit size (S-33100) provides sufficient material to



**Figure 1.** Normalized fluorescence excitation and emission spectra of SYBR Safe DNA gel stain, determined in the presence of DNA.

stain ~20 minigels; the 4 L unit size (S-33101) provides sufficient material to stain ~80 minigels. The stain is provided in 0.5X TBE buffer (45 mM Tris-borate, 1 mM EDTA, pH ~8.3). For compact storage and ease of dispensing, the 4 L unit size comes in a cube-shaped container with a removable spigot.

The SYBR Safe DNA Gel Stain Starter Kit (S-33110) is a convenient packaging of the 1 L unit size plus one SYBR Safe photographic filter (S-37100).

#### Storage

SYBR Safe DNA gel stain is stable for at least 6 months when stored at room temperature. The stain is provided in a light-proof container and should not be stored in a secondary container where it could be exposed unnecessarily to light.

#### Handling and Disposal

SYBR Safe DNA gel stain showed no or very low mutagenic activity when tested by an independent, licensed testing laboratory, and this stain is not classified as hazardous waste under U.S. Federal regulations. The safety testing included three well-established mammalian cell-based tests (Table 1), a battery of well-established Ames-test bacterial strains (Figure 2) and extensive testing for environmental safety (Table 2). Nevertheless, please exercise appropriate care and judgment when using this reagent, and dispose of the stain in compliance with all pertaining local regulations.

**Table 1.** Summary of mammalian cell-based tests for DNA genotoxicity.

Test *	Cell Type	Test Result with S9 Activation †	Test Result without S9 Activation †
Transformation test <sup>1</sup>	Syrian hamster embryo (SHE) cells	Not applicable	Negative
Chromosomal aberration test <sup>2</sup>	Cultured human peripheral blood lymphocytes	Negative	Negative
Forward-mutation test <sup>3,4</sup>	L5178Y TK mouse lymphoma cells	Negative	Negative

\* All tests were performed by Covance Laboratories, Inc., Vienna, VA, an independent testing laboratory. † S9, a mammalian extract obtained from Aroclor™ 1254-induced rat liver.


1. Fundamental and Molecular Mechanisms of Mutagenesis 356:1 (1996); 2. Evans, H.J., in *Chemical Mutagens, Principles and Methods for their Detection Vol 4*, A. Hollaender, Ed., Kluwer Academic/Plenum Publishers (1976) pp. 1–29; 3. Mutation Res 72, 447 (1980); 4. Mutation Res 59, 61 (1979).

## Protocols

### Staining Nucleic Acids after Electrophoresis (Method 1)

**1.1 Soak the gel in SYBR Safe stain.** Place the gel in a plastic container, such as a pipet-tip box lid or a household food-storage container. Do not use a glass container, as the dye in the staining solution may adsorb to the walls of the container and result in poor gel staining. Add sufficient SYBR Safe DNA gel stain to cover the gel. A 50 mL volume is sufficient for staining most standard minigels. To stain larger gels, increase the volume of staining solution in proportion to the increased gel volume, and ensure that the entire gel is fully immersed during staining.

**1.2 Incubate for 30 minutes.** Protect the gel and staining solution from light by covering it with aluminum foil or by placing it in the dark. Gently and continuously agitate the gel at room temperature (e.g., on an orbital shaker at 50 rpm). No destaining is required.



**Figure 2.** Summary of Ames test results for mutagenicity. Samples were pre-treated with a mammalian S9 fraction and then tested using the indicated Ames test strain. With strains TA97a, TA98, TA100 and TA102, a result of less than twofold-above-background suggests that the compound is nonmutagenic in the test; whereas, a result of greater than this value suggests that the compound is mutagenic in the test. With strains TA1535, TA1537 and TA1538, a result of less than threefold-above-background suggests that the compound is nonmutagenic in the test; whereas, a result of greater than this value suggests that the compound is mutagenic in the test. All tests were performed by Covance Laboratories, Inc., Vienna, VA, an independent testing laboratory.

### Precasting SYBR Safe Stain in Agarose Gels (Method 2)

**2.1 Prepare the agarose gel directly in SYBR Safe DNA gel stain.** SYBR Safe stain is provided in 0.5X TBE buffer. Simply substitute SYBR Safe stain for 0.5X TBE when preparing the molten agarose. The agarose/SYBR Safe stain mixture may be heated in the microwave. As with precasting gels with ethidium bromide, the mobility of nucleic acid fragments in the gel may be somewhat slower when run in these gels, compared to their mobility in the gel without stain.

**2.2 Run the gel normally.** No post-staining or destaining is needed.

### Viewing and Photographing the Gel

Stained gels can be viewed using a standard 300 nm transilluminator, a 254 nm epi- or transilluminator or a blue-light transilluminator, such as the Clare Chemical DarkReader™ transilluminator. DNA stained with SYBR Safe stain can also be visualized and analyzed using imaging systems equipped with an excitation source in the UV range or between 470–530 nm. Refer to the excitation/emission characteristics of SYBR Safe stain (Figure 1) in selecting the optimal filter sets to use, or contact the instrument manufacturer for advice.

Stained gels can be photographed using Polaroid® 667 black-and-white print film and SYBR Safe photographic filter (S-37100). Molecular Probes' SYPRO® photographic filter (S-6656) or a

**Table 2.** Summary of environmental safety test results.

Analysis *	Method	Results
Aquatic toxicity	Fathead minnow CA Title 22 acute screening	Not classified as hazardous or toxic to aquatic life
Ignitability	EPA 1010	Not ignitable (>212°F)
Corrosivity	EPA 150.1	Not corrosive (pH = 8.25)
Corrosivity (by Corrositex®)	DOT-E 10904	Category 2 noncorrosive
Reactivity	EPA 9010B/9030A	No reactivity detected

\* All tests were independently confirmed by AMEC Earth and Environmental San Diego Bioassay Laboratory, San Diego, CA.

Kodak® Wratten #9 filter will also work well. Using this film and one of these filters, SYBR Safe DNA gel stain provides approximately twice the detection sensitivity as ethidium bromide using a photographic filter appropriate for ethidium bromide. A

standard ethidium bromide photographic filter is not appropriate for use with SYBR Safe DNA gel stain. Gels stained with SYBR Safe stain can also be imaged using a CCD camera or a laser-based scanner.

---

**Product List** *Current prices may be obtained from our Web site or from our Customer Service Department.*

Cat #	Product Name	Unit Size
S-33100	SYBR Safe™ DNA gel stain .....	1 L
S-33101	SYBR Safe™ DNA gel stain .....	4 L
S-33110	SYBR Safe™ DNA Gel Stain Starter Kit *with 1 L of SYBR Safe™ DNA gel stain (S-33100) and one photographic filter (S-37100)* .....	1 kit
S-37100	SYBR Safe™ photographic filter.....	each

---

**Contact Information**

Further information on Molecular Probes' products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Leiden, the Netherlands. All others should contact our Technical Assistance Department in Eugene, Oregon.

Please visit our Web site — [www.probes.com](http://www.probes.com) — for the most up-to-date information.

**Molecular Probes, Inc.**

29851 Willow Creek Road, Eugene, OR 97402  
Phone: (541) 465-8300 • Fax: (541) 335-0504

**Customer Service:** 6:00 am to 4:30 pm (Pacific Time)  
Phone: (541) 335-0338 • Fax: (541) 335-0305 • [order@probes.com](mailto:order@probes.com)

**Toll-Free Ordering for USA and Canada:**  
Order Phone: (800) 438-2209 • Order Fax: (800) 438-0228

**Technical Assistance:** 8:00 am to 4:00 pm (Pacific Time)  
Phone: (541) 335-0353 • Toll-Free (800) 438-2209  
Fax: (541) 335-0238 • [tech@probes.com](mailto:tech@probes.com)

**Molecular Probes Europe BV**

Poortgebouw, Rijnsburgerweg 10  
2333 AA Leiden, The Netherlands  
Phone: +31-71-5233378 • Fax: +31-71-5233419

**Customer Service:** 9:00 to 16:30 (Central European Time)  
Phone: +31-71-5236850 • Fax: +31-71-5233419  
[eurorder@probes.nl](mailto:eurorder@probes.nl)

**Technical Assistance:** 9:00 to 16:30 (Central European Time)  
Phone: +31-71-5233431 • Fax: +31-71-5241883  
[eurotech@probes.nl](mailto:eurotech@probes.nl)

Molecular Probes' products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

Molecular Probes' SYBR Safe™ DNA gel stain is licensed under existing patents solely for use in staining nucleic acids in gels.

Several of Molecular Probes' products and product applications are covered by U.S. and foreign patents and patents pending. Our products are not available for resale or other commercial uses without a specific agreement from Molecular Probes, Inc. We welcome inquiries about licensing the use of our dyes, trademarks or technologies. Please submit inquiries by e-mail to [busdev@probes.com](mailto:busdev@probes.com). All names containing the designation ® are registered with the U.S. Patent and Trademark Office.

Copyright 2004, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.